

PrimeStop Cap

The PrimeStop Cap filter functions as a microbial barrier¹

Needle-free valve

With finger stopper and finger grip designed to reduce the risk of touch contamination.

Air Vent Filter in the Spike


Vented spike with air vent filter. The bacterial and viral filter paper efficiency is higher than 99.99%!²

- Newly designed grip plate with valves made of Tritan™, integrated back check valves and needle-free connections.
- The valve-protected connection is a closed system.³

AirStop Filter

Designed to protect against air infusion by preventing the infusion line from running dry. The filter is available in Cyto-Set® Infusomat Space and Cyto-Set® Infusion.⁴

Easy re-fill of drip-chamber due to short tubing above and low priming volume of grip-plate.



Cyto-Set® IV Administration Set

B. Braun offers a convenient solution for the application of cytostatic drugs in several situations. The new Cyto-Set design protects users in accordance with our Safe Infusion Therapy concept.

Cyto-Set® IV Administration Set

Cyto-Set®	Product	Type	Light protection	AirStop	Priming volume cm ³ = ml	Length in cm	Code no. (REF) Box of 20
Preparation							
	Cyto-Set® Mix	with 1 needle-free valve	–	–	3.3	32	A2900N
		with 1 needle-free valve, with 0.2 µm filter	–	–	5.6	32	A2903N
	Cyto-Set® Line	with 1 needle-free valve	■	–	3.3	32	A2906N
		without needle-free valve	–	–	2.8	30	A2581NF
		without needle-free valve, with 0.2 µm filter	–	–	5.3	32	A2582NF
		Application (Gravity)					
	Cyto-Set® Infusion	with 3 needle-free valves	–	–	26.5	215	A1687
		with 5 needle-free valves	–	–	26.8	215	A1686SNF
		with 3 needle-free valves	–	■	26.6	215	A1688
		with 5 needle-free valves	–	■	26.8	215	A1687SNF
Application (Pump) Infusomat® Space							
	Cyto-Set® Infusomat® Space	with 3 needle-free valves	–	–	27.3	210	8250917SP
		with 3 needle-free valves	–	■	27.3	210	835917SP
		with 5 needle-free valves	–	–	27.6	215	8250817SP
		with 5 needle-free valves	–	■	27.6	215	835817SP
		with 3 needle-free valves	■	–	27.3	210	8250920SP
		with 3 needle-free valves	■	■	27.3	210	835920SP
		with 5 needle-free valves	■	–	27.6	210	8250820SP
		with 5 needle-free valves	■	■	27.6	210	835820SP
		with 5 needle-free valves, with 0.2 µm filter	–	–	31.0	225	8250414SP
		with 5 needle-free valves, with 0.2 µm filter	–	■	31.0	225	835414SP
	Cyto-Set® Pump Adapter	4 needle-free valves	–	–	3.0	15	A1673SO

References

1. Dr. rer. nat. J. Brünke (2018). Test Report-Microbial Barrier Testing by means of Bacterial Aerosol Contamination. Quality Labs BT GmbH Nurember, Report 20018-0006.
2. Dr. Heyer (2015). Laboratory Test Report: BFE / VFE Filter paper 6739680. Data on file.
3. Abstract - Closed system test by means of Sodium Fluorescein signed by Dr. rer. nat. J. Brünke Quality Labs BT GmbH Nuremberg, Report 1816.2, 01.09.2015.
4. Pheripherenöse Schwerkraftinfusionen - Intrafix® SafeSet mit Vorteilen gegenüber herkömmlichen Infusionssystemen written by lic. rer. pol. Andreas Frei, Die Schwester Der Pfleger 43. Jahrg. 5/04

B. Braun of Canada Ltd. | Mississauga, ON
1-855-882-7286 | BBraun.ca