

Perifix[®] FX Springwound Epidural Catheters, Sets and Trays

Perifix FX Benefits


- Helps reduce paresthesia and vessel puncture rates
- Flexible
- Easy to place and remove
- Consistent
- Reliable



Easier to advance to minimize the incidence of catheter repositioning for better analgesia quality



Consistent feel and performance with minimal variation in tip deflection



Reduced paresthesia and vessel puncture with a soft and atraumatic distal tip



Easy insertion and removal with reduced threat of breakage or uncoiling with the spring embedded in polyamide



Closed-tip, multi-port design with three lateral holes for optimal drug spread, fewer unilateral blocks and a potentially higher detection rate of intravascular catheter placement

Also available in an open-tip, uni-port design

Ordering Information

PERIFIX FX Springwound Epidural Catheters

19 Ga. x 39 in. (100 cm) Radiopaque Polyamide Springwound PERIFIX FX Catheter – open tip, with catheter connector and threading assist guide (not MRI compatible)

Product Code	REF No.	Case Qty.
EC190F	333514	25

19 Ga. x 39 in. (100 cm) Radiopaque Polyamide Springwound PERIFIX FX Catheter – closed tip and three lateral side ports, with catheter connector and threading assist guide (not MRI compatible)

EC19CF	333512	25
--------	--------	----

PERIFIX FX Springwound Continuous Epidural Sets – 17 Ga. Needle

17 Ga. X 3½ in. (90 mm) Winged Tuohy Needle Set with 19 Ga. Springwound open tip catheter, catheter connector and threading assist guide

CE17TF	332234	12
--------	--------	----

17 Ga. X 3½ in. (90 mm) Winged Tuohy Needle Set with 19 Ga. Springwound closed tip catheter, catheter connector and threading assist guide

CE17TFC	332238	12
---------	--------	----

Continued on other side...

Call 1-800-227-2862 to order or visit www.perifixfx.bbraunusa.com/FX for more information.

