

TEC-1000

Tamper-Evident Syringe Cap

Your Pharmacy's Security System



B. Braun's TEC-1000, Tamper-Evident Syringe Cap, was designed with patient safety in mind. The TEC-1000 allows the healthcare worker to easily verify the integrity of the syringe by visual inspection. When the TEC-1000 is securely placed on an IV luer lock syringe and verified by the healthcare worker, unused syringes may be returned to the pharmacy to be re-dispensed. The TEC-1000 assists your staff in enhancing medication safety and reducing waste and associated costs.

B | BRAUN
SHARING EXPERTISE

TEC-1000

Tamper-Evident Syringe Cap

- Ease of Use...



CAP

Cap the filled IV luer lock syringe with a twist of the hand. Cap cannot be removed without breaking the protective barrier.



INSPECT

The clear outer shell allows visual inspection to verify medication integrity. If the shell is detached, the pharmacist, nurse, and/or patient will know that the integrity of the syringe's contents may have been compromised.



UNLOCK

The cap easily unlocks at dosage time by pulling the clear outer shell straight off and removing the red syringe cap by twisting counter-clockwise.

- Patient Safety...

The pharmacist can feel confident that the integrity of the filled syringe is protected.

- Security...

Ideal for use with syringes containing controlled substances or chemo drugs.

- Economy...

Unused medication syringes can be safely returned to the pharmacy and re-dispensed, reducing the cost of medication waste.

- Convenience...

Unique blister package of 10 protects the sterility of unopened caps.

- Quality Control...

Supports your pharmacy's efforts to provide a high quality product to the patient.

Ordering Information

Product Code	Reorder Number	Description	Case Quantity
TEC-1000	418004	Tamper-Evident Syringe Cap	500

B | BRAUN
SHARING EXPERTISE

www.bbraunusa.com

Rx only

© 2002 B. Braun Medical Inc.
IVX09-rev05/02